

Setiprant

Catalog No: #S3520


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Setiprant
Brief Description	Inhibitors
Purification	98.00%
Target Name	GPR antagonist
Calculated MW	402.42
Formulation	C24H19FN2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: C1CN(CC2=C1N(C3=C2C=C(C=C3)F)CC(=O)O)C(=O)C4=CC=CC5=CC=CC=C54 Pathways: Endocrinology/Hormones Receptor: CRTh2 Boiling pt: Melting pt: Solubility: DMSO: 36 mg/mL Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Diamant Z., et al. Setiprant, a selective CRTH2 antagonist, reduces allergen-induced airway responses in allergic asthmatics. Clin Exp Allergy. 2014 Aug;44(8):1044-52.

Note: This product is for in vitro research use only