

Halofuginone hydrobromide

Catalog No: #S3524


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Halofuginone hydrobromide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	495.59
Formulation	C16H17BrClN3O3 HB
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Others SMILES:C1C[C@H]([C@@H](NC1)CC(=O)CN2C=NC3=CC(=C(C=C3C2=O)Cl)Br)O.Br Pathways:Others Receptor:ProRS Boling
 pt:Melting pt:Solubility:DMSO Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath
 for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhou H. et al. ATP-directed capture of bioactive herbal-based medicine on human tRNA synthetase. Nature. 2013 Feb 7;494(7435):121-4.

Note: This product is for in vitro research use only