

5S rRNA modifier

Catalog No: #S3532


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	5S rRNA modifier
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	176.17
Formulation	C9H8N2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CC1=C(C(=O)N2C=CN=C2)C=CO1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 52 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Spitale RC, et al. RNA SHAPE analysis in living cells. Nat Chem Biol. 2013 Jan;9(1):18-20.

Note: This product is for in vitro research use only