

## Apilimod mesylate

Catalog No: #S3533


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Apilimod mesylate
Brief Description	Inhibitors
Purification	99.00%
Target Name	IL Receptor inhibitor
Calculated MW	610.701
Formulation	C25H34N6O8S2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology

systemSMILES:CC1=CC=CC(=C1)/C=N/NC2=NC(=NC(=C2)N3CCOCC3)OCCC4=CC=CC=N4.CS(=O)(=O)O.CS(=O)(=O)O  
 Pathways: Immunology/Inflammation  
 Receptor: IL-12; IL-23  
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM/mL  
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Wada Y, et al. Selective abrogation of Th1 response by STA-5326, a potent IL-12/IL-23 inhibitor. Blood. 2007 Feb 1;109(3):1156-64.

Note: This product is for in vitro research use only