

## TETRABENAZINE

Catalog No: #S3534


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	TETRABENAZINE
Brief Description	Inhibitors
Purification	98.00%
Target Name	Monoamine Transporter inhibitor
Calculated MW	317.429
Formulation	C19H27NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
 SMILES: CC(C)CC1CN2CCC3=CC(=C(C=C3C2CC1=O)OC)OC  
 Pathways: Membrane transporter/Ion channel  
 Receptor: VMAT-2  
 Boiling pt: Melting pt: Solubility: DMSO : 58 mg/mL  
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Jankovic, J. and J. Beach, Long-term effects of tetrabenazine in hyperkinetic movement disorders. *Neurology*, 1997. 48(2): p. 358-62.

Note: This product is for in vitro research use only