

## SHP099 hydrochloride

Catalog No: #S3544


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	SHP099 hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	phosphatase inhibitor
Calculated MW	388.72
Formulation	C16H20Cl3N5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Others SMILES: CC1CC(=CCC2CC(CC3(O2)CC(C(O3)C(=CC(C)C)C)C)OP(=O)(O)OC(=O)C4C=C(C(C5C4(C(=CC=C1)CO5)O)O)C)C Pathways: Others Receptor: SHP-2 Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Yellow powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Chen YN, et al. Nature. 2016 Jul 7;535(7610):148-52.

Note: This product is for in vitro research use only