


## Description

Product Name	LX-4211
Brief Description	Inhibitors
Purification	99.00%
Target Name	SGLT inhibitor
Calculated MW	424.93
Formulation	C21H25ClO5S
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:MetabolismSMILES:CCOC1=CC=C(C=C1)CC2=C(C=CC(=C2)[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)SC)O)O)O)ClPathways:GPCR/G  
ProteinReceptor:SGLT1; SGLT2Boiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:White powderRemark:For obtaining a higher solubility ,  
please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Zambrowicz B, et al. Clin Pharmacol Ther. 2012 Aug;92(2):158-69.;Powell DR, et al. J Pharmacol Exp Ther. 2013 May;345(2):250-9.;

Note: This product is for in vitro research use only