

Epacadostat

Catalog No: #S3548


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Epacadostat
Brief Description	Inhibitors
Purification	100.00%
Target Name	IDO inhibitor
Calculated MW	436.99
Formulation	C ₁₁ H ₁₃ BrFN ₇ O ₄ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: C1=CC(=C(C=C1N/C(=C\2/C(=NON2)NCCNS(=O)(=O)N)/N=O)Br)F
 Pathways: Metabolism
 Receptor: IDO1
 Boiling pt: Melting pt: Solubility: Soluble in DMSO
 Appearance: White powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
 Stock solution can be stored below -20 °C for several month.

References

1. Vacchelli E, et al. Oncoimmunology. 2014 Dec 15;3(10):e957994.

Note: This product is for in vitro research use only