

Lvguidingan

Catalog No: #S3563


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lvguidingan
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	272.1735
Formulation	C13H15Cl2NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: CCC(C)NC(=O)/C=C/C1=CC(=C(C=C1)Cl)Cl
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting
 Solubility: Soluble in DMSO
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
 Stock solution can be stored below -20 °C for several months.

References

1. Patejdl R, et al. Phenytoin inhibits contractions of rat gastrointestinal and portal vein smooth muscle by inhibiting calcium entry. *Neurogastroenterol Motil.* 2015 Oct;27(10):1453-65.

Note: This product is for in vitro research use only