


## Description

Product Name	SU 5205
Brief Description	Inhibitors
Purification	99.00%
Target Name	VEGFR inhibitor
Calculated MW	239.2494
Formulation	C15H10FNO
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:O=C1NC2=C(C=CC=C2)/C1=C/C3=CC=C(F)C=C3 Pathways:Angiogenesis; Tyrosine

Kinase/AdaptorsReceptor:VEGFR2Boiling pt:Melting pt:Solubility:Soluble in DMSO Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Tian XF, et al. PACE4 regulates apoptosis in human pancreatic cancer Panc 1 cells via the mitochondrial signaling pathway. Mol Med Rep. 2016 Dec;14(6):5205-5210.

Note: This product is for in vitro research use only