


Description

Product Name	NSC 247030
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	290.143
Formulation	C15H9Cl2NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Others SMILES: C1=CC=C2C(=C1)/C(=C/C3=CC(=C(C=C3)Cl)Cl)/C(=O)N2 Pathways:Others Receptor:Others Boiling pt:Melting pt:Solubility:Soluble in DMSO Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only