

Pyridoxal phosphate


Catalog No: #S3578

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pyridoxal phosphate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others antagonist
Calculated MW	247.14
Formulation	C8H10NO6P
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Drug-induced dyskinesia
 SMILES: CC1=NC=C(C(=C1O)C=O)COP(=O)(O)O
 Pathways: Others
 Receptor: Metabolism
 Boiling pt: Melting pt: Solubility: Water: soluble
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Modak MJ. Pyridoxal 5' phosphate: a selective inhibitor of oncornaviral DNA polymerases. *Biochem Biophys Res Commun.* 1976 Jul 12;71(1):180-7.

Note: This product is for in vitro research use only