

Src I 1


Catalog No: #S3593

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Src I 1
Brief Description	Inhibitors
Purification	99.00%
Target Name	Src inhibitor
Calculated MW	373.4
Formulation	C22H19N3O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:COC1=C(C=C2C(=C1)C(=NC=N2)NC3=CC=C(C=C3)OC4=CC=CC=C4)OCPathways:AngiogenesisReceptor:Lck;
SrcBoiling pt:Melting pt:Solubility:DMSO: 6.2 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Tian G, et al. Structural determinants for potent, selective dual site inhibition of human pp60c-src by 4-anilinoquinazolines. *Biochemistry*. 2001 Jun 19;40(24):7084-91.

Note: This product is for in vitro research use only