

## Cutamesine dihydrochloride

Catalog No: #S3597


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cutamesine dihydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Sigma receptor agonist
Calculated MW	441.43
Formulation	C <sub>23</sub> H <sub>32</sub> N <sub>2</sub> O <sub>2</sub> ·2HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: COC1=C(C=C(C=C1)CCN2CCN(CC2)CCCC3=CC=CC=C3)OC.Cl.Cl Pathways: GPCR/G

Protein Receptor:  $\sigma$ 1 Boiling pt: Melting pt: Solubility: Water: 100 mM; DMSO: 20 Mm Appearance: White powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Tagashira H, et al. Biochim Biophys Acta. 2013, 1830(4):3082-94.

Note: This product is for in vitro research use only