

Benodanil


Catalog No: #S3601

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|------------------------------------------------|
| Product Name | Benodanil |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Others |
| Calculated MW | 323.13 |
| Formulation | C13H10INO |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:OthersSMILES:C1=CC=C(C=C1)NC(=O)C2=CC=CC=C2IPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Chloroform, DMSOAppearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only