

D-Pantothenic acid sodium salt

Catalog No: #S3602


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	D-Pantothenic acid sodium salt
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	241.22
Formulation	C ₉ H ₁₆ NNaO ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CC(C)(CO)[C@H](C(=O)NCCC(=O)[O-])O.[Na+] Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in water Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only