

Gastrodenol

Catalog No: #S3605


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Gastrodenol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	704.47
Formulation	C ₁₂ H ₁₀ BiK ₃ O ₁₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C(C(=O)O)C(CC(=O)[O-])(C(=O)O)O.C(C(=O)O)C(CC(=O)[O-])(C(=O)O)[O-].[K].[K].[K].[Bi+3] Pathways: Others Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 40 mg/mL Appearance: White powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

N/A

Note: This product is for in vitro research use only