

Gamma-Oryzanol

Catalog No: #S3606


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Gamma-Oryzanol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	602.9
Formulation	C40H58O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular system/ Digestive system/

ImmunologySMILES:CC(CCC=C(C)C)[C@H]1CC[C@@]2([C@@]1(CCC34C2CCC5[C@]3(C4)CC[C@@H](C5(C)C)OC(=O)/C=C/C6=CC(=C(C=C6)O)OC)C)CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in ChloroformAppearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Rao YP, et al. Biochem Biophys Res Commun. 2016 Oct 28;479(4):747-752.

Note: This product is for in vitro research use only