

Kitasamycin

Catalog No: #S3607


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Kitasamycin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	701.85
Formulation	C35H59NO13
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C[C@@H]1C/C=C/C=C/[C@@H]([C@@H](C[C@@H]([C@@H]([C@H]([C@@H](CC(=O)O1)O)OC)O[C@H]2[C@@H]([C@@H]([C@@H]([C@H](O2)C)O[C@H]3C[C@@]([C@H]([C@@H](O3)C)O)(C)O)N(C)C)O)CC=O)O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in ethanol Appearance: White powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhang C, et al. Anaerobe. 2016 Dec;42:188-196.

Note: This product is for in vitro research use only