

Dioxopromethazine

Catalog No: #S3609


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dioxopromethazine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	316.42
Formulation	C17H20N2O2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CC(CN(C)CC1=CC=CC=C1S(=O)(=O)C2=CC=CC=C21)N(C)C Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only