

## Sodium Glucuronic Acid

Catalog No: #S3611


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sodium Glucuronic Acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	147.15
Formulation	C9H6FN
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:[Na+].O[C@H]1[C@H](O)[C@@H](OC(O)[C@@H]1O)C([O-])=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in waterAppearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only