

Carbocisteine

Catalog No: #S3612


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Carbocisteine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	179.19
Formulation	C5H9NO4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:C(C(C(=O)O)N)SCC(=O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Slight soluble in WaterAppearance:White powderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Takeda K, et al. Br J Pharmacol. 2010 Jul; 160(6): 1399C1407.

Note: This product is for in vitro research use only