

## Ivosidenib

Catalog No: #S3617


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ivosidenib
Brief Description	Inhibitors
Purification	99.48%
Target Name	Dehydrogenase inhibitor
Calculated MW	582.97
Formulation	C28H22ClF3N6O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Cancer SMILES: C1CC(=O)N([C@@H]1C(=O)N(C2=CC(=CN=C2)F)[C@@H](C3=CC=CC=C3Cl)C(=O)NC4CC(C4)(F)F)C5=NC=CC(=C5)C#N  
 Pathways: Metabolism Receptor: IDH1 Boiling pt: Melting pt: Solubility: DMSO: 39 mg/mL Appearance: White powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Agresta Samuel V, et al. From PCT Int. Appl. (2015), WO 2015127172 A1 20150827.

Note: This product is for in vitro research use only