

Solcitinib

Catalog No: #S3620


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Solcitinib
Brief Description	Inhibitors
Purification	99.82%
Target Name	JAK inhibitor
Calculated MW	389.18
Formulation	C22H23N5O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: O=C(C1CC1)NC2=NN3C(C4=CC=C(C(N5CC(C)(C)C5)=O)C=C4)=CC=CC3=N2
 Pathways: Angiogenesis; Chromatin/Epigenetic; JAK/STAT signaling; Stem Cells
 Receptor: JAK1
 Boling pt: Melting pt: Solubility: DMSO: 20 mg/ML
 Appearance: White powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kahl L, et al. Lupus. 2016 Apr 6. pii: 0961203316640910.;

Note: This product is for in vitro research use only