

## CP-640186 hydrochloride

Catalog No: #S3622


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	CP-640186 hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Acetyl-CoA Carboxylase inhibitor
Calculated MW	522.08
Formulation	C30H35ClN3O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Others SMILES: C1C[C@H](CN(C1)C2CCN(CC2)C(=O)C3=C4C=CC=CC4=CC5=CC=CC=C53)C(=O)N6CCOCC6. CIP Pathways: Proteases/Proteasome Receptor: ACC1; ACC2 Boiling pt: Melting pt: Solubility: DMSO: 48 mg/ml Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Harwood HJ Jr, et al. J Biol Chem. 2003 Sep 26;278(39):37099-111.;

Note: This product is for in vitro research use only