


Description

Product Name	IQ-1S
Brief Description	Inhibitors
Purification	97.92%
Target Name	JNK inhibitor
Calculated MW	247.25
Formulation	C15H9N3O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: O=N=C1C2=C(C3=NC4=CC=CC=C4N=C31)C=CC=C2 Pathways: MAPK Receptor: JNK1; JNK2; JNK3 Boiling pt: Melting pt: Solubility: DMSO: 32 mg/mL Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Igor A. Schepetkin et al. J Pharmacol Exp Ther, 2015 Jun, 353(3): 505-516.

Note: This product is for in vitro research use only