

Madrasin

Catalog No: #S3628


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Madrasin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	311.34
Formulation	C16H17N5O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: O=C1NC(NC2=NC(C)=C3C=CC(OC)=CC3=N2)=NC(C)=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Andrea Pawellek et al. J Biol Chem, 2014 Dec 12, 289(50): 34683-34698.

Note: This product is for in vitro research use only