

AZD-9291 mesylate?

Catalog No: #S3634


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	AZD-9291 mesylate?
Brief Description	Inhibitors
Purification	99.46%
Target Name	EGFR inhibitor
Calculated MW	595.71
Formulation	C29H37N7O5S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:CN1C=C(C2=CC=CC=C21)C3=NC(=NC=C3)NC4=C(C=C(C(=C4)NC(=O)C=C)N(C)CCN(C)C)OC.CS(=O)(=O)OPathways:Angiogenesis; JAK/STAT signaling; Tyrosine Kinase/AdaptorsReceptor:EGFR (L858R); EGFR (L858R/T790M)Boiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Patent WO 2013014448 A1: 2 - (2, 4, 5 - substituted -anilino) pyrimidine derivatives as egfr modulators useful for treating Y .

Note: This product is for in vitro research use only