

Pifithrin-β hydrobromide

Catalog No: #S3637


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pifithrin-β hydrobromide
Brief Description	Inhibitors
Purification	99.24%
Target Name	p53 inhibitor
Calculated MW	349.29
Formulation	C ₁₆ H ₁₇ BrN ₂ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CC1=CC=C(C=C1)C2=CN3C4=C(CCCC4)SC3=N2.Br Pathways: Apoptosis Receptor: p53 Boiling pt: Melting pt: Solubility: DMSO: 26 mg/mL Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Abdelalim EM, et al. Biochem Biophys Res Commun. 2012 Apr 13;420(3):605-10.

Note: This product is for in vitro research use only