

## Leukadherin-1

Catalog No: #S3638


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Leukadherin-1
Brief Description	Inhibitors
Purification	98.00%
Target Name	Integrin activator
Calculated MW	421.49
Formulation	C22H15NO4S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:C1=CC=C(C=C1)CN2C(=O)C(=C/C3=CC=C(O3)C4=CC=C(C=C4)C(=O)O)/SC2=SPathways:Cytoskeletal SignalingReceptor:CD11b/CD18Boiling pt:Melting pt:Solubility:DMSO: 6 mg/mL (Ultrasonic)Appearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Reed JH, et al. J Biol Chem. 2013 Mar 29;288(13):9077-83.

Note: This product is for in vitro research use only