

Endoxifen E-isomer hydrochloride

Catalog No: #S3645


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Endoxifen E-isomer hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Estrogen/progestogen Receptor inhibitor
Calculated MW	409.95
Formulation	C25H28ClNO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: CC/C(=C(/C1=CC=C(C=C1)O)\C2=CC=C(C=C2)OCCNC)/C3=CC=CC=C3.Cl Pathways: Endocrinology/Hormones Receptor: ERR

Boling pt: Melting pt: Solubility: DMSO: 30 mg/mL Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Chae YJ, et al. Eur J Pharmacol. 2015 Apr 5;752:1-7.

Note: This product is for in vitro research use only