


Description

Product Name	MX69
Brief Description	Inhibitors
Purification	97.27%
Target Name	Mdm2 inhibitor
Calculated MW	474.57
Formulation	C27H26N2O4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES: Cc1cc(ccc1C)NS(=O)(=O)c2ccc3NC(C4CC=CC4c3c2)c5ccc(cc5)C(=O)O
 Pathways:ApoptosisReceptor:Mdm2Boiling pt:Melting pt:Solubility:DMSO: 95 mg/mL (200.18 mM); Ethanol: 41 mg/mL (86.39 mM)Appearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Gu L, et al. Y Cell. 2016, 30(4):623-636.

Note: This product is for in vitro research use only