


Description

Product Name	APY29
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	332.36
Formulation	C17H16N8
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:C1CC1C2=CC(=NN2)NC3=NC(=NC=C3)NC4=CC5=C(C=C4)N=CN5Pathways:OthersReceptor:IRE1αBoiling pt:Melting pt:Solubility:DMSO: 10 mg/mLAppearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wang L, et al. Nat Chem Biol. 2012 Dec;8(12):982-9.;

Note: This product is for in vitro research use only