

## Eliglustat

Catalog No: #S3663


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

|                   |  |
|-------------------|--|
| Product Name      | Eliglustat                                     |
| Brief Description | Inhibitors                                     |
| Purification      | 99.00%   |
| Target Name       | Transferase inhibitor                          |
| Calculated MW     | 404.54   |
| Formulation       | C23H36N2O4                                     |
| Storage           | 3 years -20°C powder;2 years -80°C in solvent; |

## Images



## Product Description

## Research

Area:MetabolismSMILES:CCCCCCCC(=O)N[C@H](CN1CCCC1)[C@@H](C2=CC3=C(C=C2)OCCO3)OPathways:MetabolismReceptor:glucosylceramide synthaseBoiling pt:Melting pt:Solubility:DMSO: 80 mg/mL (197.75 mM);Ethanol: 80 mg/mL (197.75 mM)Appearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Shayman JA. Drugs Future. 2010 Aug 1;35(8):613-620.

Note: This product is for in vitro research use only