

Tectoridin?

Catalog No: #S3665

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tectoridin?
Brief Description	Inhibitors
Purification	99.00%
Target Name	Reductase inhibitor
Calculated MW	462.4
Formulation	C22H22O11
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research

Area:OthersSMILES:COC1=C(C=C2C(=C1O)C(=O)C(=CO2)C3=CC=C(C=C3)O)O[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)OPathway s:Endocrinology/Hormones; EnzymeReceptor:Aldose reductaseBoiling pt:Melting pt:Solubility:Pyridine, Methanol, EthanolAppearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Li Y, et al. Luminescence. 2016 Feb;31(1):38-46.

Note: This product is for in vitro research use only