

## Atractylodin

Catalog No: #S3666


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Atractylodin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Lipid inhibitor
Calculated MW	182.22
Formulation	C13H10O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:C/C=C/C#CC#C/C=C/C1=CC=CO1

Pathways:MetabolismReceptor:lipaseBoiling pt:Melting pt:Solubility:Chloroform, Dichloromethane, Ethyl Acetate, DMSO,

AcetoneAppearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Tao J, et al. Zhongguo Zhong Yao Za Zhi. 2011 Jan;36(1):27-31.

Note: This product is for in vitro research use only