

## Dehydroandrographolide Succinate Potassium Salt

Catalog No: #S3667


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Dehydroandrographolide Succinate Potassium Salt
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	570.67
Formulation	C28H35OK10
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C[C@@]12CC[C@H]([C@@]1[C@H]1CCC(=C)[C@H]2/C=C/C3=CCOC3=O)(C)COC(=O)CCC(=O)[O-]OC(=O)CCC(=O)O.[K+]

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Hsieh MJ, et al. Oncotarget. 2015 Oct 13;6(31):30831-49.

Note: This product is for in vitro research use only