

Galangin

Catalog No: #S3668


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Galangin
Brief Description	Inhibitors
Purification	99.00%
Target Name	ERK inhibitor; NF-κB inhibitor
Calculated MW	270.2
Formulation	C15H10O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others
SMILES: C1=CC=C(C=C1)C2=C(C(=O)C3=C(C=C(C=C3O)O)O)O
Pathways: MAPK; NF-κB Receptor: ERK; NF-κB
Boiling pt: Melting pt: Solubility: Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone
Appearance: Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Wang Y, et al. Toxicology. 2014 Dec 4;326:9-17.

Note: This product is for in vitro research use only