

## Curcumenol


Catalog No: #S3669

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Curcumenol
Brief Description	Inhibitors
Purification	98.00%
Target Name	P450 inhibitor
Calculated MW	234.3
Formulation	C15H22O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: CC1=C(C)C2=C(C)C(O)C3=C(C)C=C(C)C3O2 Pathways: Metabolism Receptor: CYP450 Boiling pt: Melting pt: Solubility: Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Hamdi OA, et al. Int J Mol Sci. 2015 Mar 6;16(3):5180-93.

Note: This product is for in vitro research use only