

Mollugin

Catalog No: #S3673


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mollugin
Brief Description	Inhibitors
Purification	98.00%
Target Name	HER inhibitor; JAK inhibitor
Calculated MW	284.3
Formulation	C17H16O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CC1(C=CC2=C(O1)C3=CC=CC=C3C(=C2C(=O)OC)O)C Pathways: Angiogenesis; Chromatin/Epigenetic; JAK/STAT signaling; Tyrosine Kinase/Adaptors; Stem Cells Receptor: HER2; JAK2 Boiling pt: Melting pt: Solubility: Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Zhang L, et al. Biochem Biophys Res Commun. 2014 Jul 18;450(1):247-54.

Note: This product is for in vitro research use only