

Hydroxysafflor yellow A

Catalog No: #S3674


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hydroxysafflor yellow A
Brief Description	Inhibitors
Purification	97.00%
Target Name	TGF-beta/Smad inhibitor
Calculated MW	612.53338
Formulation	C27H32O16
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: C1=CC(=CC=C1/C=C/C(=C2/C(=C(C(=O)C(C2=O))([C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)O)[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)O)O)O Pathways: Stem Cells Receptor: TGFβ Boling pt: Melting pt: Solubility: Pyridine, Methanol, Ethanol Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhu HJ, et al. Cytotechnology. 2015 Oct;67(5):885-92.

Note: This product is for in vitro research use only