

Sesamol?

Catalog No: #S3676


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sesamol?
Brief Description	Inhibitors
Purification	99.90%
Target Name	AMPK activator
Calculated MW	138.12
Formulation	C7H6O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:C1OC2=C(O1)C=C(C=C2)OPathways:PI3K/Akt/mTOR signalingReceptor:AMPKBoling pt:Melting pt:Solubility:Chloroform, Dichloromethane, Ethyl Acetate, DMSO,?Appearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wu XL, et al. Inflamm Res. 2015 Aug;64(8):577-88.

Note: This product is for in vitro research use only