

## Tebanicline

Catalog No: #S3694


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Tebanicline
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR agonist
Calculated MW	198.65?
Formulation	C9H11ClN2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
SMILES: C1CN[C@H]1COC2=CN=C(C=C2)Cl  
Pathways: Neuroscience  
Receptor: nAChR  
Boiling pt: Melting pt:  
Solubility: Soluble in DMSO  
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.  
Stock solution can be stored below -20 °C for several month.

## References

1. Hayashi T, et al. *Neurosci Lett.* 2016 Dec 7;638:46-50.

Note: This product is for in vitro research use only