

Pepstatin

Catalog No: #S3695


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pepstatin
Brief Description	Inhibitors
Purification	99.09%
Target Name	Amino Acids and Derivatives inhibitor
Calculated MW	685.89
Formulation	C34H63N5O9
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:C[C@@H](C(=O)N[C@@H](CC(C)C)[C@H](CC(=O)O)O)NC(=O)C[C@@H]([C@H](CC(C)C)NC(=O)[C@H](C(C)C)NC(=O)[C@H](C(C)C)NC(=O)CC(C)C)OPathways:OthersReceptor:aspartic proteaseBoiling pt:Melting pt:Solubility:DMSO: 100 mg/mL (145.79 mM)Appearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Seelmeier S, et al. Proc Natl Acad Sci U S A. 1988, 85(18), 6612-6616.

Note: This product is for in vitro research use only