

polygalaxanthone III

Catalog No: #S3696


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	polygalaxanthone III
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	568.39
Formulation	C25H28O15
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:COC1=C(C=C2C(=C1)C(=O)C3=C(C(=C(C=C3O2)O)[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO[C@H]5[C@@H]([C@](

CO5)(CO)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol,

EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lv C, et al. J Sep Sci. 2014 May;37(9-10):1103-10.

Note: This product is for in vitro research use only