

## Lasmiditan

Catalog No: #S3704


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Lasmiditan
Brief Description	Inhibitors
Purification	99.00%
Target Name	5-HT Receptor agonist
Calculated MW	377.14
Formulation	C19H18F3N3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
 SMILES: CN1CCC(CC1)C(=O)C2=NC(=CC=C2)NC(=O)C3=C(C=C(C=C3F)F)F  
 Pathways: GPCR/G Protein;  
 Neuroscience Receptor: 5-HT1F  
 Boiling pt: Melting pt: Solubility: Soluble in DMSO  
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Nelson DL, et al. Cephalalgia. 2010 Oct;30(10):1159-69.

Note: This product is for in vitro research use only