

Mivebresib

Catalog No: #S3712


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mivebresib
Brief Description	Inhibitors
Purification	99.00%
Target Name	Epigenetic Reader Domain inhibitor
Calculated MW	459.11
Formulation	C22H19F2N3O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer
 SMILES: CCS(=O)(=O)NC1=CC(=C(C=C1)OC2=C(C=C(C=C2)F)F)C3=CN(C(=O)C4=C3C=CN4)C
 Pathways: Chromatin/Epigenetic Receptor: BET
 Boiling pt: Melting pt: Solubility: DMSO: 37 mg/mL
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. EJ Faivre et al. Y Research, 2016

Note: This product is for in vitro research use only