

# Rolapitant


Catalog No: #S3716

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                   |  |
|-------------------|--|
| Product Name      | Rolapitant   |
| Brief Description | Inhibitors   |
| Purification      | 98.00%   |
| Target Name       | Neurokinin 1 antagonist  |
| Calculated MW     | 500.48   |
| Formulation       | C <sub>25</sub> H <sub>26</sub> F <sub>6</sub> N <sub>2</sub> O <sub>2</sub> |
| Storage           | 3 years -20°C powder; 2 years -80°C in solvent;                              |

## Images



## Product Description

### Research

Area: Others SMILES: C[C@H](C1=CC(=CC(=C1)C(F)(F)F)C(F)(F)F)OC[C@]2(CC[C@]3(CCC(=O)N3)CN2)C4=CC=CC=C4 Pathways: Neuroscience Receptor: NK1 Boiling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Rapoport B, et al. Eur J Y. 2016 Apr;57:23-30.

Note: This product is for in vitro research use only