

Daminozide

Catalog No: #S3719


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Daminozide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Histone Demethylase inhibitor
Calculated MW	160.17
Formulation	C6H12N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC(CCC(NN(C)C)=O)=O
Pathways: Chromatin/Epigenetic Receptor: KDM2/7
JmJCB
Boiling pt: Melting pt: Solubility: Soluble in DMSO
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Reed DJ, Moore TC, Anderson JD. Science. 1965 Jun 11;148(3676):1469-71.

Note: This product is for in vitro research use only