

Tebanicline hydrochloride

Catalog No: #S3720


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tebanicline hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR inhibitor
Calculated MW	235.11
Formulation	C ₉ H ₁₂ Cl ₂ N ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: C1CN[C@H]1COC2=CN=C(C=C2)Cl
 Pathways: Neuroscience
 Receptor: nAChR
 Boiling pt: Melting pt: Solubility: Soluble in DMSO
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Nirogi R, Jabar SL, Jayarajan P, et al. Eur J Pharmacol. 2011 Oct 1;668(1-2):155-62.

Note: This product is for in vitro research use only